



MEDICAL INSTRUMENTS

Vision

Actively improve the quality of services and products to win customer satisfaction, while for each employee to bring material and spiritual aspects of the harvest, so that real wealth growth, fulfill the social responsibility and benefit society.

Mission

To fill each corner of the world with love.
To possess health for each family worldwide.

Core Values

Honesty and Faith: Honesty and Faithful character and conduct.

Optimism: A positive and optimistic attitude.

Responsibility: The spiritual sense of responsibility.

Integrity: Integrity and selfless mind.

Tolerance: Broad-mindedness of grace.

Universal love: Spirit of loving the whole world.

Philosophy

Integrity first, service supreme, unlimited love.

Management Concept

Institutionalization, standardization, humanism and internationalization.

Enterprise Concept

Three "emotions" : passion, enthusiasm and sincere.
Four "heart" : a sense of responsibility, gratitude heart, self-confidence, subject to heart.

Service Concept

Respect for customers, understanding customers, continued to provide quality products and services, to become lifelong partner of our clients.

Talent Concept

Attitude first, no excuses, no impossible, integrity, responsibility.

Work Concept

Work hard, go all out, sincere communication, courtesy, focus on learning, fun work.

Conviction

Keep their promises, highly responsible, practice what they preach, lead by example, says it will do.

Social Values

With love to do business, with action back to society, to benefit more people.

Team Concept

United, selfless, love and support.

Contents

SPHYGMOMANOMETER
Page 01~08

STETHOSCOPE
Page 09~20

THERMOMETER
Page 21~30

PHYSICIAN SCAN
Page 31~34

MOBILITY AID
Page 35~42

HOSPITAL EQUIPMENT
Page 43~56





SPHYGMOMANOMETER

1. material of valve: Copper Alloy or Zinc Alloy
2. material of tube: PVC or Latex, ID:3.00mm, OD:8.00mm, Length48-50cm
3. Gauge: Zinc Alloy, cassette mechanism: copper or stainless steel, weight:130-135g
4. material of bladder: PVC or latex,
infant size: 11X6cm, child size: 18X8.5cm, adult size: 22X12cm
5. material of cuff:nylong or cotton,
infant size: 28X8cm, child size: 34.5X11cm, adult size: 50X14cm
6. Maximum error tolerance +/- 3 mm Hg
7. Any colors of tube and cuff is available

ANOMETER
RY
RIGHT
ARM



GlenMedicon LLC

Add: 7122 S Sheridan Rd Ste 2 - 1012Tulsa, OK 74133-2748
 Tel: 1-918-8091339
 Fax: 1-918-523-7121
 E-mail: info@glenmedicon.com
 Http://www.glenmedicon.com

Aneroid Sphygmomanometer

- 300mmHg Gauge
- Cotton or Nylon cuff
- With or without D-ring is available
- Artery label
- Vinyl or Nylon zippered carrying case



GM001-101

Aneroid Sphygmomanometer



GM001-110

Aneroid Sphygmomanometer With Separate Single Head Stethoscope



GM001-111

Aneroid Sphygmomanometer with Attached Single Head Stethoscope



GM001-120

Palm Type Aneroid Sphygmomanometer Plastic Shell 65mm Diameter dial



GM001-121

Palm Type Aneroid Sphygmomanometer Metal Shell 65mm Diameter dial

*Best service to all old and new clients!
 Best greeting to all patients in the world!*



GM001-122
Palm Type Aneroid Sphygmomanometer
Metal Shell
50mm Diameter dial



GM001-130
Arm Type
Aneroid Sphygmomanometer



GM001-140
White plastic cover Type
Aneroid Sphygmomanometer



GM001-123
Palm Type Aneroid Sphygmomanometer
Metal Shell
65mm Diameter dial



GM001-124
Palm Type Aneroid Sphygmomanometer
Metal Shell
50mm Diameter dial



GM001-125
Palm Type Aneroid Sphygmomanometer
Plastic Shell
65mm Diameter dial



GM001-150
Desk Type
Aneroid Sphygmomanometer



GM001-160
Wall Type
Aneroid Sphygmomanometer



GM001-126
Palm Type Aneroid Sphygmomanometer
Plastic Shell
55mm Diameter dial

*Best service to all old and new clients!
Best greeting to all patients in the world!*

GM001-127
Five Size Blood Pressure Kit

1PC Palm Type Aneroid Sphygmomanometer Gauge
5Pcs Sizes Cuff with latex baladder & quick release
luer adapter:

- 1PC Infant Cuff: 28cm x8cm
- 1PC Child Cuff: 34.5cm x11cm
- 1PC Adult Cuff: 50.8cm x14cm
- 1PC Large Cuff: 63cm x17.5cm
- 1PC Thigh Cuff: 72cm x21cm
- 1Pc Penlight
- 1PC Dual head Stethoscope
- 1Pc Neurological Hammer
- 1Pc Nylon Carrying Case

Form of packing:
10pcs/carton
Size of packing:
Export carton:600x420x300mm




GM001-170
Floor Type Aneroid Sphygmomanometer

Form of packing:
One machine fill in a white box, and every one boxes in one export carton.

Size of packing:
White box:495x460x170mm
Export carton:500x470x185mm

Mercury Sphygmomanometer

300mmHg scale glass tubes
Metal case
Cotton or nylon cuff,
With or without D-ring is available
Artery label
Latex bladder with double tubes



GM001-200
Professional Mercury Sphygmomanometer
Desk type



GM001-20
1 LCD Display



GM001-20
2 Light Display



GM001-210
Mercury Sphygmomanometer
Wall Type



GM001-220
Mercury Sphygmomanometer
Floor Type

*Best service to all old and new clients!
Best greeting to all patients in the world!*

Accessories



GM001-401
Zipper Case



GM001-403
Bladder



GM001-404
Bulb



GM001-402
Cuff



GM001-405
Air-release Valve



GM001-406
End Valve



GM001-407 Cuff
Material: Nylon or cotton

Fully automatic electronic blood pressure monitor

Display system: Digital display system/LCD
 Measuring range: 2.7kpa-40kpa(20mmHg-300mmHg)
 Accuracy: 0.5kPa(4mmHg)
 Power source: Type AA batteries*4 for arm type
 Type AA batteries*2 for wrist type
 Operating environment: -5°C less than 95RH
 Automatic power off: To be automatically cut off after few minutes of non use to save energy



GM001-301
Arm type



GM001-302
Arm type



GM001-303
Arm type



GM001-311
Wrist type



GM001-312
Wrist type



GM001-313
Wrist type