



GlenMedicon LLC Medical Instruments Co.,Ltd.

Best Greeting To All Patients In The World !

▶ Recommend new disposable medical products

www.glenmedicon.com



GlenMedicon LLC

Add: 7122 S Sheridan Rd Ste 2 - 1012Tulsa, OK 74133-2748

Tel: 1-918-8091339

Fax: 1-918-523-7121

E-mail: info@glenmedicon.com

[Http://www.glenmedicon.com](http://www.glenmedicon.com)





www.glenmedicon.com



GlenMedicon Medical Equipment & Supplies is a minority and women owned business enterprise. We provide practitioners, patients and the general public with the highest quality products at the lowest available prices. We envision a world in which medical products can be accessible to all.

Regardless of the market or product category, our diversified relationships with all CE approved and FDA approved product manufacturing partners around the globe, can help your company drive revenue and overall valuation. Our strong track record of proven performance is the direct result of our decision to build the company's future on a specialized approach to marketing - one that promotes a unique understanding of the needs of the medical community. Our global accessibility gives us resources to effectively fulfill our commitment to excellence, now and in the years to come.

GlenMedicon successfully represents some of the most well respected Speciality Medical Device manufactures in the world today, allowing us to expand our product offering to multiple departments in the Hospitals, Diagnostic Laboratories , Rehabilitation Centers, Home health & Assisted living Centers as well as to Surgery Centers and a variety of Specialist Clinics. We specialize in respiratory, urology, wound care , catheter , syringes and needles. Our primary concern is your satisfaction.

Best Greetings to all the patients in the world.

We pray no disaster ,no flu, We hope everyone is healthy and fortunate.



GM037-400
Auto Self -destruction Syringe

Size : 0.05ml, 0.1ml, 0.2ml, 0.5ml, 1ml



GM037-500
Disposable destruction Syringe

Size: 2ml, 3ml, 5ml ,10ml





GM037-600
Safety Clip Syringe

Size: 1ml, 3ml, 5ml, 10ml, 20ml

Operation procedure

- 1 Fix the clip before drawing.
- 2 Same operation as normal syringe in drawing and injecting
- 3 After injection, push the clip over the needle and lock it onto the clip base



GM037-800
Safety Auto Destruct Insulin Syringe

Size: 100ul



GM038-900
Insulin needle

Needle Size	Needle length	Package
30G	8mm	100pcs/box
31G	8mm	100pcs/box
31G	6mm	100pcs/box
31G	5mm	100pcs/box
32G	4mm	100pcs/box
33G	4mm	100pcs/box



GM037-700
Sliding -sheath Safety Syringe

Size: 1ml, 3ml, 5ml, 10ml, 20ml



GM039-300
Combi stopper -Red

For Single Use Only
Non Toxic & Pyrogen free



GM039-400
Disposable Bite Block

GM034-103 I.V Infusion Set

Transparent & flexible big drip chamber
Available in DEHP free model
With Y site
Soft and kink resistant PVC tubing ensures uniform flow rate
Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates
Standard tube length : 150 cm



GM034-104 I.V Infusion Set with Filter

I.V. Infusion set with 0.2 micron IV in-line filter in order to prevent any particles, entrapped air and bacteria & fungi passing into the venous system
Available in DEHP free model
Transparent & flexible drip chamber
With Y site
Soft and kink resistant PVC tubing ensures uniform flow rate
Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates
Standard tube length : 150 cms



GM034-105 I.V Infusion Set with Flow Regulator

Sharp and easy piercing air vented spike
Flow regulator integrated into IV infusion set for precise flow control
Soft and kink resistant PVC tubing ensures uniform flow rate
Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates
Standard tube length : 150 cms



GM034-106 I.V Flow Regulator

Precise flow rate control of IV fluids with range of 5 to 250 ml/hr
With "Y" injection port (latex or latex free) for intermittent medication
Also available with thumb support "Y" injection port and needle free "Y" site
Male luer connector at one end and female luer connector at other end
Standard tube length : 20 cm without "Y" site & 40 cm with "Y" site



GM034-300

Disposable Infusion Pump

- The product is used in transfusion for the treatment of postoperation pain. It is composed of solution storage device, power plant, liquid filter unit, current-limiting device, accessories, option switch, etc. It can be divided into automatic pump, PCA pump and multi-functional pump.
- Stock solution: 100ml,200ml,275ml
- Flow rate: 2ml/h,3ml/h,4ml/h,5ml/h,6ml/h,8ml/h,10ml/h



GM034-300A
Multi-function



GM034-300B
PCA Self-control fluid feeding



GM034-300C
Automatic

Anesthesia Kit



GM038-800E/S



GM038-800S

Accessories	800-E/S	800-E	800-S
Epidural needle	16G/18G 80mm/90mm	Same	
Spinal needle	25G/27G / 113/123mm(Pencilpoint)		20G/22G
Epidural catheter	Pull-resistance length > 90cm	Same	
Epidural Catheter adapter	One	Same	
Liquid medicine filter	(blue)5 μ	Same	Same
Air filter	(Yellow)0.5μ	Same	
LOR Syringe	One Glass 5ml	Same	Same
Plastic Syringe	Two(3ml and 20ml)	1 pc 20ml	1 pc 5ml
negative pressure tube	one	one	
fenesrated sheet	60cm x 60cm hole 10cm	Same	Same
Small Nonwoven sheet	1 pc 50cm x 50cm	Same	Same
Big Nonwoven sheet	1pc 70cm x95cm	Same	Same
Catheter adhesive silce	One	Same	
Adhesive traps	One 5cm x 40cm , 2pc-1.25cm x 10cm	Same	
Compresses	Two 8 plyers 8cm x 10cm, one 5 x6 with cut	Same	Same
Sponge brush	Three(curving/Straight)	Same	Same
Surgical latex gloves	One pair 7.5/7.0/6.0	Same	Same
Syringe needle	One OD 0.5mm L20mm One OD 0.7mm L32mm One OD 1.6mm L38mm	Same	Same

Anesthesia Kit (Mini)



Bowie-Dick Test Pack




Accessories	900-E/S	900-E	900-S
Epidural needle	16G/18G 80mm/90mm	Same	
Spinal needle	25G/27G 113/123MM(Pencilpoint)		20G/22G/25G
Epidural catheter	Pull-resistance length > 90cm	Same	
Epidural Catheter adapter	One	Same	
Liquid medicine filter	(blue) 5 μ	Same	Same
Air filter	(Yellow) 0.5μ	Same	
LOR Syringe	One Glass 5ml	Same	Same
Plastic Syringe	Two(3ml and 20ml)	1 pc 20ml	1 pc 5ml
Syringe needle	One OD 0.5mm L20mm One OD 0.7mm L32mm One OD 1.6mm L38mm	Same	Same

Test pack is a 125mm×115mm block about 18mm thick, over wrapped with blue paper.

The interior of the pack contains blotter type paper, a single plastic foam sheet(4mm thick)&a single paper sheet with indicator ink.

- 1 Point out the presence of air leaks in the autoclave immediately and accurately;
- 2 Operation Condition: 134°C, 3.5minutes
- 3 Color Change: Yellow to Brown



GM024-200 Hemodialysis Blood Tubing Set

Hemodialysis Blood Tubing Set consists of venous lines, used during dialysis and are attached with A.V Fistula needle and dialyzer.

Features:

Medical grade materials displaying good biocompatibility Fluid Pathway is smooth inside, no clotting or air bubbles
Tubes are flexible, kink -resistant
Color coding
Pump segment has excellent anti-fatigue performance
Luer-lock connectors according to international standard
Various series, good applicability

Type	OD of Pump segment	Quality of chamber	To Be Used With
I-A	12.0mm	1	Fresenius; Gambro; Baxter;nipro;etc
I-B	9.8mm	1	B. Braun
II-A	12.0mm	2	Fresenius; Gambro; Baxter;nipro;b.braun;nikkiso;etc
II-B	9.8mm	2	Baxter
C	9.8mm	---	Fresenius; Gambro; Baxter; Nipro;etc



GM024-300 Dilation and drainage Set

Dilation and drainage set is a surgical implement which used for dilation of a tract and drainage in percutaneous nephrostomy.

Feature

The dilator is made from medical grade polymer material.
The drainage tube is made from medical grade silicone.
Good biocompatibility can minimize harm to human body.
Peeling sheath has suitable hardness and flexibility which can offer best maintenance during surgery.
Round tip of dilator can minimize trauma to human body.
Provide wide selections based on client' s requirement.

Instruction for use

- 1 By introducing a guide needle into pyelo in percutaneous nephrostomy.
 - 2 Withdrawing the core of guide needle.
 - 3 Introducing a guide wire through guide needle.
 - 4 Withdrawing the body of guide needle.
 - 5 Progressing from the smallest to the largest appropriate size by passing the dilator over the guide wire while maintaining the position the guide wire.
 - 6 The largest size dilator must be coated with peeling sheath.
 - 7 Withdrawing the guide wire and dilator, and maintaining the peeling sheath in the body. The surgery can be manipulated through peeling sheath.
 - 8 After surgery, the drainage tube can be positioned into body for drainage.
 - 9 Withdrawing peeling sheath.
- Size:Fr6- Fr30



GM024-400
Bladder puncture fistula

Feature

Big aperture can be beneficial for drainage during surgery
Good structure is good for easy control during surgery
Size: Fr.15, Fr.18, Fr.24



GM024-500
Suction Plus Endotracheal Tube
(Endotracheal Tube with Evacuation Lumen)

Suction Plus Endotracheal Tube is designed with the evacuation lumen in order to provide protection against risk of aspiration and decrease the rate of ventilation-associated pneumonia(VAP) in the hospital.
Size ID: 6.0mm,6.5mm,7.0mm,7.5mm,8.0mm,8.5mm,9.0mm



GM024-600 Guide wire



GM024-700
Endobronchial Blocker Tube

Endobronchial blocker tube is a device which can be inserted down a tracheal tube after tracheal intubation so as to block off the right or left main bronchus of the lungs in order to be able to achieve a controlled one lung ventilation.
The Lung tissue distal to the obstruction will collapse, thus allowing the surgeon' s view and access to relevant structures within the thoracic cavity.
Our blocker tubes are used to achieve lung separation and one lung ventilation as an alternative to double-lumen tubes(DLT) and serve as a choice for pediatric patients as well.
Standard PU Cuff Type
Size: 5FR,7FR,9FR



GM015-300
Tracheostomy Tube With/With Cuff

Made from non -toxic PVC
Be used in tracheostomy
90°Curvature
With 15mm standard connector
X-ray opaque line throughout the length of the tube
Available with of standard low pressure cuff and without cuff
Pre--connected neck -tape
Size: 6.0mm-10.0mm

Disposable Electrosurgical Pencil

- 1 170mm-ABS Pencil appearance
- 2 7cm Stainless steel tip
- 3 2.36mm-diameter of the tip
- 4 3m-PVC medical cable
- 5 Maximal out put current:1A
- 6 Frequency confine:0-1.0MHZ
- 7 Ways of aseptic technique:Sterile EO
- 8 Current resistance:200



GM040-500



GM040-600



GM040-200A



GM040-200B



GM040-200C



GM040-300

Disposable Medical Razor

Carbon steel or stainless steel, single or double blade.
All razors come with the highest quality stainless steel blades or carbon steel, comfort-coated for smooth shaves.
Packing: 100pcs/box, 1000pcs/carton,63x32x21cm



GM102-150
Hydrocolloid Dressing with thin border

Size : 5cm x 5cm, 7.5cm x7.5cm,
10cm x 10cm, 15cm x 15cm
15x20cm, 20cm x 20cm



GM102-160
Hydrocolloid Dressing extra thin 0.6mm

Size : 5cm x 5cm, 7.5cm x7.5cm, 10cm x 10cm,
15cm x 15cm;20cm x20cm
15x20cm, 20cm x 20cm



GM102-170
Foam Dressing Laminated with PU film

Size:5cm x 5cm, 10cm x 10cm, 10cm x20cm, 15cm x15cm, 20cm x 20cm



GM102-180
Tranparent Fixing Dressing

Size: 6cmx7cm ,10cm x12cm

GM102-190
I.V Cannula Tranparent Fixing Dressing

Size: 6cmx7cm ,10cm x12cm



GM102-200
I.V Cannula Fixing Dressing with Pad

Size: 6cmx7cm



GM102-210
Non woven Fixing Dressing with Pad

Size:8cm x 6cm ,or other sizes per customer' s demand.



GM102-220
Wound Skin Closure

Size: 1/8" x 3" (3mm x 75mm); 1/4"



GM102-230
Lattice Wound tape





GM102-240
Sterile Soluble Hemostatic Gauze

Hemostatic gauze is an all-natural material produced from plant cellulose. When used properly, hemostatic gauze stops bleeding quickly and thoroughly in a wide variety of applications.
 Size: 5cm x 7.5cm



GM102-250
Absorbent hemostatic gauze



GM096-200
Strip sports adhesive

The strip sports adhesive tape uses cotton cloth serving as a base material and adhesive cement mainly made of zinc oxide and rubber, and the adhesive cement is coated in a shape of strip. When the product is used, firstly clean and dry the skin, tear out the product according to the needed length, and apply it to the needed part. The product has the sports protection effect, but it can't be directly applied to the wounds.

Size: 5 / 7.5 / 10 / 15 / 20cm X 5M




GM096-300
Orthopedic Wool Padding

Material: 100% Cotton or Synthetic fiber
 Size: 5 / 7.5 / 10 / 15 / 20cm X 4M



GM096-400
Transparent film Rolls

PU adhesive dressing
 The PU adhesive dressing is mainly formed by two parts: a base material (PU) and separate paper. The PU base material is soft, thin and light and mainly used for fixing. The product is used for the protection of various open injuries. The product is low in sensitization, excellent in biological compatibility, bacteria-free and high in air permeation, and doesn't irritate or damage the skin nor affect the natural functions of the skin.



Material

- ① PET backing
 - ② PU Film
 - ③ Medical Pressure Sensitive adhesive
 - ④ Release paper
- Size: 5cm x 10m, 10cm x 10m, 15cm x 10m

Rapid Testing Series



Product	Format	Specimen	Packing
HCG	Strip(2.0/2.5/3.0mm)	Urine	100T/BOX
HCG	Device	Urine	25T/BOX
HCG	Midstream	Urine	20T/BOX
LH	Strip	Urine	100T/BOX
LH	Device	Urine	25T/BOX
LH	Midstream	Urine	20T/BOX
HIV 1/2	Strip	Serum/Plasma	100T/BOX
HIV 1/2	Device	Serum/Plasma	25T/BOX
HIV 1/2	Device	WB/Serum/Plasma	25T/BOX
Syphilis	Strip	Serum/Plasma	100T/BOX
Syphilis	Device	Serum/Plasma	25T/BOX
Syphilis	Device	WB/Serum/Plasma	25T/BOX
HAV IgM	Device	Serum/Plasma	25T/BOX
HBsAG	Strip	Serum/Plasma	100T/BOX
HBsAG	Device	Serum/Plasma	25T/BOX
HBsAG	Device	WB/Serum/Plasma	25T/BOX
HCV	Strip	Serum/Plasma	100T/BOX
HCV	Device	Serum/Plasma	25T/BOX
HCV	Device	WB/Serum/Plasma	25T/BOX
Typhoid IgG/IgM	Device	Serum/Plasma	25T/BOX
TOXO IgG/IgM	Device	Serum/Plasma	25T/BOX
Malaria pf/pv	Device	Whole Blood	25T/BOX
Malaria pf/pan	Device	Whole Blood	25T/BOX
Malaria pf	Device	Whole Blood	25T/BOX
Flu A&B	Strip	Nasopharyngeal swab	20T/BOX
Dengue(IgG/IgM)	Device	Whole Blood	25T/BOX
TB	Strip	Serum/Plasma	100T/BOX
TB	Device	Serum/Plasma	25T/BOX
TB	Device	WB/Serum/Plasma	25T/BOX
M.P(IgM)	Device	Serum/Plasma	25T/BOX
H.P Antigen	Device	Feces	25T/BOX
H.P Antibody	Device	Serum/Plasma	25T/BOX
Chlamydia	Device	Swab/Urine	20T/BOX
AFP	Strip	Serum/Plasma	100T/BOX
AFP	Device	Serum/Plasma	25T/BOX
FOB	Strip	Feces	100T/BOX
FOB	Device	Feces	25T/BOX
PSA	Strip	Serum/Plasma	100T/BOX
PSA	Device	Serum/Plasma	25T/BOX
TroponinI	Device	WB/Serum/Plasma	25T/BOX

Urinalysis Reagent Strips



Name	Specifications	Units	Packaging specifications	Application Analyzer	Test items
Urinalysis Reagent Strips	1 item	bottle	300 bottles / Box	visual	glucose
Urinalysis Reagent Strips	2 items	bottle	300 bottles / Box	visual	protein, glucose, ketone, pH, protein, pH, glucose
Urinalysis Reagent Strips	3 items	bottle	300 bottles / Box	visual	pH, glucose, protein
Urinalysis Reagent Strips	4 items	bottle	300 bottles / Box	visual	pH, glucose, protein, occult blood
Urinalysis Reagent Strips	8A	bottle	40 bottles / Box	MA-4210, URITEST-200	urobilinogen, bilirubin, ketone, occult blood, protein, glucose, pH, nitrite
Urinalysis Reagent Strips	8N	bottle	40 bottles / Box	visual	urobilinogen, occult blood, bilirubin, ketone, glucose, protein, pH, nitrite
Urinalysis Reagent Strips	10A	bottle	40 bottles / Box	CLINITEK ANALYZER	leukocytes, nitrite, urobilinogen, protein, pH, occult blood, specific gravity, ketone, bilirubin, glucose
Urinalysis Reagent Strips	10U	bottle	40 bottles / Box	URITEST-100/200/300	leukocytes, ketone, urobilinogen, bilirubin, specific gravity, occult blood, protein, glucose, nitrite, Ph
Urinalysis Reagent Strips	11U	bottle	40 bottles / Box	URITEST-100/200/300	pH, nitrite, glucose, ascorbic acid, specific gravity, occult blood, protein, bilirubin, urobilinogen, ketone, leukocytes
Urinalysis Reagent Strips	11H	bottle	40 bottles / Box	HY-100	urobilinogen, occult blood, bilirubin, ketone, leukocytes, glucose, protein, pH, nitrite, specific gravity, ascorbic acid

Name	specifications	Units	Packaging specifications	Application Analyzer	Test items
Urinalysis Reagent Strips	11D	bottle	40 bottles / Box	MAC-120, YX-200, etc	leukocytes, nitrite, urobilinogen, protein, pH, occult blood, specific gravity, ascorbic acid, ketone, bilirubin, glucose
Urinalysis Reagent Strips	11T	bottle	40 bottles / Box	MA-4280	ascorbic acid, leukocytes, occult blood, urobilinogen, bilirubin, pH, specific gravity, ketone, protein, nitrite, glucose
Urinalysis Reagent Strips	10B	bottle	40 bottles / Box	BM-200	glucose, protein, urobilinogen, bilirubin, ketone, nitrite, occult blood, leukocytes, pH, specific gravity
Urinalysis Reagent Strips	11A	bottle	40 bottles / Box	HY-100, HL3460B	leukocytes, urobilinogen, occult blood, bilirubin, ketone, glucose, protein, pH, nitrite, specific gravity, ascorbic acid
Urine Analyzer	HY-100	piece	1 piece/box		



GM109-120

GM109-110

Saline Nasal Rinse Kit



GM159-108

Vaginal Speculum with Light



GM159-110

Vomit/Emesis bag

Features: a one way valve that prevents leaking or spilling even if the bag is dropped; includes printed graduations for measurement.

GM159-110: PE +Paper

GM159-120: PP+ PE material ,1000CC, Packing: 4 sleeves of 25 Each



GM159-120



GM159-130 NEEDLE COUNTER

Description
NEEDLE COUNTER 20CT FM BL ADH
NEEDLE COUNTER 20CT DBL MAG
NEEDLE COUNTER 60CT DBL MAG
NEEDLE COUNTER 30CT DBL MAG
NEEDLE COUNTER 30CT FM BL MAG
NEEDLE COUNTER 20CT FM BL MAG
NEEDLE COUNTER 40CT FM BL ADH
NEEDLE COUNTER 10CT MAG ADH
NEEDLE COUNTER 60CT DBL MAG DEEP
NEEDLE COUNTER 30CT DBL FM ST
NEEDLE COUNTER 80CT DBL FM
NEEDLE COUNTER 60CT MAG ADH
NEEDLE COUNTER 60CT DBL FM ST
NEEDLE COUNTER 20CT FM BL SM BOX
NEEDLE COUNTER 15CT MAG CLR LID

Description
NEEDLE COUNTER, 10-COUNT, FOAM BLOCK WITH ADHESIVE
NEEDLE COUNTER, 10-COUNT, FOAM BLOCK WITH MAGNET
NEEDLE COUNTER, 10-COUNT, MAGNET, CLEAR TOP
NEEDLE COUNTER, 15- COUNT, FOAM STRIPS WITH MAGNET
NEEDLE COUNTER, 15-COUNT, MAGNET, CLEAR TOP
NEEDLE COUNTER, 20-COUNT, FOAM BLOCK WITH MAGNET
NEEDLE COUNTER, 20-COUNT, DOUBLE MAGNET
NEEDLE COUNTER, 40-COUNT, FOAM BLOCK WITH ADHESIVE
NEEDLE COUNTER, 40-COUNT, FOAM BLOCK WITH MAGNET
NEEDLE COUNTER With /Without blade remover, 60-COUNT, DOUBLE MAGNET