



**GlenMedicon LLC**

Add: 7122 S Sheridan Rd Ste 2 - 1012Tulsa, OK 74133-2748 Tel:  
1-918-8091339  
Fax: 1-918-523-7121  
E-mail: [info@glenmedicon.com](mailto:info@glenmedicon.com)  
[Http://www.glenmedicon.com](http://www.glenmedicon.com)

# GlenMedicon LLC

DISPOSABLE  
HOSPITAL UNIFORM SERIES

2019





**GlenMedicon** Medical Equipment & Supplies is a minority and women owned business enterprise. We provide practitioners, patients and the general public with the highest quality products at the lowest available prices. We envision a world in which medical products can be accessible to all.

Regardless of the market or product category, our diversified relationships with all CE approved and FDA approved product manufacturing partners around the globe, can help your company drive revenue and overall valuation. Our strong track record of proven performance is the direct result of our decision to build the company's future on a specialized approach to marketing - one that promotes a unique understanding of the needs of the medical community. Our global accessibility gives us resources to effectively fulfill our commitment to excellence, now and in the years to come.

**GlenMedicon** successfully represents some of the most well respected Speciality Medical Device manufactures in the world today, allowing us to expand our product offering to multiple departments in the Hospitals, Diagnostic Laboratories , Rehabilitation Centers, Home health & Assisted living Centers as well as to Surgery Centers and a variety of Specialist Clinics. We specialize in respiratory, urology, wound care , catheter , syringes and needles.

Our primary concern is your satisfaction.

Best Greetings to all the patients in the world.

We pray no disaster ,no flu, We hope everyone is healthy and fortunate.





**GM058-101 Round Nurse Cap** 一次性圆帽

**Description:**  
Nurse Caps Manufactured of Spunbonded Polypropylene Lightweight, Cool, Comfortable Material – Soft Elastic Headband Suited To Lab, industry, Hospital And Daily Home Use.

**Specification:**  
Non Woven PP Material :  
10g,12g; 14g;16g;18g;  
Available In 18" ,21" or 24" Sizes,  
White, Green or Blue, Pink, Yellow, Red... etc.  
Sterile Method : None or Eo Sterile  
Packing: 100pcs/Poly Bag 10Bags/Carton  
carton size: 240x240x340mm



Type A: Single Elastic

**GM058-201 Clip Cap** 一次性折叠帽

**Description:**  
Manufactured From High Quality Non-woven Polypropylene High Quality Material For Increased Durability Convenient Dispenser Lightweight, Cool, Comfortable Material – Soft Elastic Headband ,and Comfortable Fold Type Suited To Lab, Industry, Hospital and Daily Home Use.

**Specification:**  
Non Woven PP Material:10g;12g;14g;16g, etc  
Color: White, Green, Blue , Pink ,Red, Yellow  
Size:20",21" ,24"  
Type: Single Elastic /Double Elastic  
Sterile Method : None or EO Sterile  
Packing: 100pcs/Poly Bag 10Bags/Carton  
carton size: 530x240x350mm



Type B: Double Elastic



Type A: Ties

**GM058-301 Doctor Cap** 一次性医生帽

**Description:**  
Manufactured From High Quality, Soft, Non-woven PP Material, High Quality Material For Increased Durability Convenient Dispenser Lightweight, Cool, Comfortable Material – Soft Elastic Headband Suited To Lab, Industry, Hospital and Daily Home Use.

**Specification:**  
Non-woven PP Material; 20g,25g,30g,40g  
Color: Green ,Blue  
Type: Tie or Elastic Back  
Size: 61x14cm ; 64x23cm etc.  
Latex Free  
The SMS ,Spun Lace Fabric is available  
Packing: 100pcs/ Poly Bag, 10Bags/Carton



Type B: Elastic Back

**GM058-401 Net Cap** 一次性网眼帽

**Specification:**  
Nylon Net Fabric  
Color: Blue  
Size: 18",21" ,24"  
Packing: 100Pcs/ Box ,  
10Box/Carton

**Size of packing:**  
Inner box:280x260x100mm  
Export carton : 510x290x530mm  
Gross. Weight: 6.0 Kgs  
Net. Weight: 5.0 Kgs



**GM058-501 Surgical Hood** 一次性太空帽

**Description:**  
A Disposable Surgical Hood Which Conforms To The Wearer's Head, and is Usable Either with a Traditional Surgical Mask Alone or With Such a Mask In Combination With a Laminar Air Flow Face Plate, Contains Bacterial Fallout From The Wearer's Head and Helps Maintain Sterile Conditions in Hospital Operating Rooms.

**Specification:**  
Non Woven PP Material: 16g,20g,25g,30g,etc  
Size: 30x124cm  
Color: White, Green, Blue  
Packing: 100Pcs/Poly Bag ; 10Bags/Carton



**GM050-601 Procedure Cap** 一次性工作帽

**Description:**  
Manufactured from High Quality, soft, Non-Woven PP material, High Quality Material for Increased Durability Convenient Dispenser Lightweight, cool, comfortable material – soft elastic headband Suited to Lab, industry, hospital and daily home use.

**Specification:**  
Non-woven PP material;20g,25g,30g,40g  
Color: Green ,Blue Type: elastic back  
Size:61x14cm ;64x23cm etc. Latex free  
The SMS ,Spunlace fabric is available  
**Form of packing:** 100pcs/poly bag, 10bags/carton  
**Size of packing:**Export carton: 50x50x30cm or 18x30x30cm  
G.W: 7.0kgs N.W: 6.0kgs



**GM058-701 SBP Head cover** 一次性套头帽

Manufactured from High Quality Non-Woven Polypropylene High Quality Material for Increased Durability Convenient Dispenser Lightweight, cool, comfortable material ; Suited to Lab, industry, hospital and daily home use.

**Specification:**  
Non woven PP material:10G;12G;14G;16G etc  
Color: White, green, blue , pink ,red, yellow  
Form of packing: 100pcs/poly bag 10bags/carton  
carton size: 450x280x260mm  
G.W.: 4.0 KGS; N.W.: 3.5 KGS



**GM058-801 Snood cap** 一次性发网帽

**Description:**  
Manufactured from High Quality, soft, Non-Woven PP material, High Quality Material for Increased Durability Convenient Dispenser Lightweight, cool, comfortable material – soft elastic headband Suited to Lab, industry, hospital and daily home use.

**Specification:**  
Non-woven PP material;10gsm/ 12gsm/13gsm/14gsm/16gsm,etc.  
Color: white, green, blue,yellow, red, etc.  
Standard size: 62x9 cm or other size per customer'demand.  
Form of packing: 100pcs/ poly bag,10bags/carton



**一次性无纺口罩**

Conforms to American NELSON N95 test standards.  
Mini-fiber high electrostatic filter, filtration efficiency of dust & powder is over 95%.  
Repel continually contamination of dangerous blood and saliva continually.  
Inner cover is well sweat-absorbed, while outer cover is water resistant.  
polypropylene cloth, low breathing resistance and high filtration efficiency.

**GM059-101 Non-woven Face Mask**

**Specification:**

Ear-Loop Type,  
Ply: 1ply, 2-Ply or 3 Ply  
Detailed Weight: 1 Ply:38g/m<sup>2</sup>, 2Ply:25g/m<sup>2</sup>+18g/m<sup>2</sup>;  
3Ply:25g/m<sup>2</sup>+20g/m<sup>2</sup>+16g/m<sup>2</sup>;  
Size: 18x9cm  
Color: White, Green,Blue, Red, Pink, Yellow, etc.  
Packing: 50Pcs/Box 40Boxes/Carton



**GM059-102 Non-woven Face Mask**

**Specification:**

With ties type  
Ply: 2-ply or 3 ply  
Detailed weight: 2ply:25G/M2; 3ply:25G/M2+20G/M2+16G/M2;  
Size: 18x9cm  
Color:white,Green,Blue,Red, Pink, Yellow, etc.  
Packing: 50pcs/box, 40 boxes/carton,  
Carton size: 520x380x380CM

**GM059-103 Non-woven Face Mask with Charcoal**

**Description:** 一次性无纺口罩 (活性炭)

Conforms to American NELSON N95 test standards.  
3-ply polypropylene cloth, low breathing resistance and high filtration efficiency.  
Adjustable nose seal, easy to fold and carry. Mini-fiber high electrostatic filter,  
filtration efficiency of dust & powder is over 95%.  
Repel continually contamination of dangerous blood and saliva continually.  
Suggested applications: school, restaurant, surgery, electronic, theaters,  
beauty salon...etc public places.

**Specification:**

Ear-loop type: 3 ply  
Detailed weight: 3ply:25G/M2+20G/M2+16G/M2;  
Size: 18x9cm Color:white,green,blue etc.  
Packing: 50pcs/box, 40 boxes/carton, Carton size: 520x380x300CM




**GM059-200 Nurse Paper Face Mask 一次性纸口罩**

**Description:**

Our Most Economical Face Masks, High Quality Wet-Strength Paper  
Comfortable Elastic Ear Loops For Easy Donning Product Description

**Specification:**


2 Ply or 1 Ply  
Size: 20x7cm  
Color: White  
Suggested Applications: Polluted Environment,  
Especially Good For Absorbing and Filtrating The Exhausted Gas  
From Motor Vehicles.  
Packing: 100 Pcs/Box 100Boxes/Carton



**GM059 -104 4-ply Anti-penetration Nonwoven Face Mask**

**四层防渗透无纺布口罩**

Size: 9 2(18)x18cm  
Material: 18g/m<sup>2</sup> Wet-laid nonwovens (or 22g/m<sup>2</sup> PP)+anti-penetration  
Punching membrane + 24g/m<sup>2</sup> melting-blown PP + 25g/m<sup>2</sup> calendering PP  
Packing: 50pcs/box, 6box/carton or 20box/carton  
Function: 4-ply structure, effectively prevent from bacteria and dust, high  
efficiency in resisting penetration, BFE up to 99%  
Production technique: ultrasonic machine-sewed, with tie and nosepiece  
Application: suitable for casualty wards, infection wards, surgery rooms,  
medical sanitation, scientific experiment, chemical industry and pharmacy



**GM059 -121 Face Mask 三层海绵层口罩**

Material: 18g/m<sup>2</sup> we-laid nonwovens (or 22g/m<sup>2</sup> PP)+ 24g/m<sup>2</sup> melting-blown PP + 25g/m<sup>2</sup> calendering PP  
Packing: 50pcs/box, 6box/carton or 20box/carton  
Function: 3-ply structure, effectively prevent from bacteria and dust, BFE up to 99%  
Production technique: ultrasonic machine-sewed, with tie and nosepiece  
Application: suitable for casualty wards, infection wards, surgery rooms, medical sanitation, scientific experiment, chemical industry and pharmacy

**GM059 -120 4-ply Non-penetrated Anti-fog Face Mask**

**四层防渗透防雾片口罩**

Size: 9.2 (18) x 18 cm, with anti-fog piece  
Material: 18g/m<sup>2</sup> wet-laid nonwovens (or 22g/m<sup>2</sup> PP) + non-  
penetrated punching membrane + 24g/m<sup>2</sup> melting-blown PP + 25g/m<sup>2</sup>  
calendering PP  
Packing: 25pcs/box, 10box/carton  
Function: 4-ply structure, effectively prevent from bacteria and dust,  
protect eyes and face from fluid, provide better safe property  
Production technique: ultrasonic machine-sewed, with tie and anti-  
fog nosepiece  
Application: suitable for casualty wards, infection wards, and scientific  
experiment



**GM059-120 镜片口罩**

**DUST-PROOF FACE MASK 一次性防尘口罩**

**Terylene Series 涤纶系列**

Low breathing resistance with high filtration efficiency. Adopted the Terylene material, which is without any toxic chemicals or stimulant. Adjustable contoured aluminum nosepiece and a soft inner foam nose cushion help to ensure proper fit and increase worker comfort. Special 3-ply design respirator produced from high filtrable and electrostatic material for resisting normal toxic and tiny micro particle. Suggested applications: chemurgy, spray-paint/printing, dyeing, garage, factories of electronic aid and food, disinfecting and the insalubrious workshops.



GM059-301A



GM059-301B  
With Respritor Valve



GM059-302A



GM059-302B  
With Respritor Valve



GM059-303A



GM059-303B  
With Respritor Valve

**DUST-PROOF FACE MASK 一次性防尘口罩**

**Compound Series 合成系列**

Low breathing resistance with high filtration efficiency. Adopted the polypropylene material, which is without any toxic chemicals or stimulant. Adjustable contoured aluminum nosepiece and a soft inner foam nose cushion help to ensure proper fit and increase worker comfort. Special 3-ply design respirator produced from high filtrable and electrostatic material for resisting normal toxic and tiny micro particle. Suggested applications: chemurgy, spray-paint/printing, dyeing, garage, factories of electronic aid and food, disinfecting and the insalubrious workshops.



GM059-305A



GM059-305B  
With Respritor Valve



GM059-305D  
With Respritor Valve

**GM059-304 N95 Face Mask (with Valve) N95 口罩**

**Specification:**

Conforms To International Standarrds:  
EN149:2001 FFP1 FFP2 CE0194  
AS/NZS 1716: 1994 NIOSH N95  
Low Breathing Resistance With High Filtration Efficiency. Adopted The Polypropylene Material, Which is Without Any Toxic Chemicals or Stimulant. Adjustable Contoured Aluminum Nosepiece and a Soft Inner Foam Nose Cushion Help To Ensure Proper Fit And Increase Worker Comfort. Easy-Breathing Exhalation Valve With Inner Silicon-Piece Can Facilitate The Air Exchange and Increase Worker Comfort. Special 3-Ply Design Respirator Produced From High Filtrable and Electrostatic Material For Resisting Normal Toxic and Tiny Micro Particle.

**Suggested Applications:**

Chemurgy, Spray-Paint/Printing, Dyeing, Garage, Factories of Electronic Aid and Food, Disinfecting and The Insalubrious Workshops.

Form of Packing: 20Pcs/ Inner Box ; 20Boxes/Carton



Type A



Type B  
With Respritor Valve



**GM059-300 MELT-BLOWN SERIES**

Conforms to FFP1 & FFP2 带铝鼻梁条系列  
Fiberglass free, odorless, non-irritating to the skin  
Adjustable contoured aluminum nosepiece helps to ensure proper fit.  
Low cost and easy to use  
Suggested applications: factory, garden, process, maintain, house cleaning, motorcycle riders, wood, and any toxic free protection.



**GM059-140 Beard Mask (Cover)**

Specification: 一次性胡须套  
Color: White, Blue, Green, etc  
Material: Non woven PP fabric  
Gram: 10G/12G/14G/16G/18G/20G etc.  
Standard Size 48 x 24CM



GM059-306



GM059-307


**FOLDING SERIES** 新型折叠防尘口罩系列

Conforms to FFP1 & FFP2 EN149: 2001 CE0194 standard.  
Folding designing, wide open for easy breathing and fold easily to fit in pocket. Special 3-ply design respirator produced from high filtrable and electrostatic material for resisting normal toxic and tiny micro particle. Adjustable contoured aluminum nosepiece permits individual facial adjustment. Suggested applications: mining, building, grinding drugs, casting, agriculture, farming, electronic, food processing, cement textile, disintegration, carving, heavy metal pollution, gardening, glass...etc.



**GM060-100 PE Apron** 一次性 PE 围裙

Specification:  
Made of Polyethylene  
With Adjustable Tie  
Available In Different Size and Color Thickness  
0.018mm, 0.020mm, 0.025mm Size : 110x70cm;  
135x71cm; 160x90cm, etc. Color: White or Blue  
Package : 100Pcs/Bag, 1000Pcs/Carton Packing: 10Pcs/Poly Bag 1000Pcs/Carton



**GM059-130 Face Guard**

Specification:  
Conforms To American Nelson N95 Test Standards.  
Ear-Loop Type: 3-Ply or 2Ply Non Woven,  
Polypropylene Cloth,  
Low Breathing Resistance  
and High Filtration Efficiency,  
and PVC Transparent Shield  
Mini-Fiber High Electrostatic Filter,  
Filtration Efficiency of Dust & Powder is over 95%.  
Repel Continually Contamination of  
Dangerous Blood and Saliva Continually.  
Inner Cover is Well Sweat-Absorbed,  
While Outer Cover is Water Resistant.  
Packing: 25Pcs/Inner Box; 20 Boxes/Carton

GM059-130  
镜片口罩



**GM061-100 PE Oversleeve** 一次性 PE 袖套

Specification:  
Size 20x40cm, 42x22cm, 46x22cm  
Thickness 0.02 mm / 0.018mm  
Made of PE. With Comfortable Elastic, Water Repellen,  
Elasticized at Both Ends to Prevent Slippage  
Fluid Resistant  
Effective Protection For Wrist and Forearm  
Clear, White, Blue, Red Colors etc  
Packing: 100Pcs/Poly Bag 2000Pcs/Carton



**GM061-101 Non woven Oversleeve** 一次性无纺袖套

Specification:  
Size 20x40cm, 42x22cm, 46x22cm  
Weight: 30/35/40 gsm at clinets' option.  
Material: Non woven PP,  
Elasticized at both ends to prevent slippage  
Fluid resistant  
Effective protection for wrist and forearm  
White, Blue, Green, Red colors etc  
Form of Packing: 100pcs/poly bag 1000pcs/carton



**GM062-101 PE Shoe Cover** 一次性 PE 鞋套

**Specification:**

Color: Blue or Green  
Material: PE Fabric: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g,  
Made of PE. With Comfortable Elalstic Opening.  
Available in Different Thickness and Color  
Size: 36x15cm, 41x15cm, 40x16cm  
Thickness: 0.03mm, 0.04mm, 0.18mm at Clients' option  
Packing: 100Pcs/Poly Bag 5000Pcs/Carton



**GM062-102 CPE Shoe Cover** 一次性 CPE 鞋套

**Specification:**

Color: Blue or Green .  
Material: CPE Fabric: 1.8g, 2.0g, 2.1g, 2.5g, 3.0g  
Made of CPE, Embossed Surface .  
with Comfortable Elastic Opening.  
Available In Different Thickness and Color  
Size: 36x15cm, 41x15cm, 40x16cm  
thickness: 0.02mm, 0.025mm, 0.28mm, 0.035mm  
At Clients' Option Packing: 100Pcs/Poly Bag 5000Pcs/Carton



**GM062-103 Non-woven Shoe Cover** 一次性无纺鞋套

**Specification:**

Color: White, Blue  
Material: PP Non-woven  
**TYPE A: PP Regular Shoe Covers** 不防滑  
Comfortable Elastic Hem, Helps You Keep Your Cool.  
Easily Worn Under A Sterile Gown. One-size-fits-most Sizing  
Reduces Your Inventory  
**TYPE B: PP Non-skid Shoe Covers** 防滑  
Non-skid Shoe Covers are Designed For Those Situations  
Where Wet or Slippery Conditions Exist  
A Proprietary Plastic Resin Is Applied To The Bottom of  
The Cover to Improve The Wearers Traciton On The Floor Weight:  
30/40/45gsm at Clients' option  
Size: 36x15cm, 38x18cm, 40x16cm  
Packing: 100Pcs/Bag 1000Pcs/Carton



**GM062-104A Boot cover** 一次性靴套

**Specification:**

Color: white, blue, green, etc.  
Material: Non-woven PP  
Comfortable elastic hem, helps you keep your cool.  
Easily worn under a sterile gown.  
Weight: 30G/35G/40G/45G.. etc.  
Standard size: 40 x35cm, or other size per customer' demand Packing:  
100pcs/poly bag, 1000 pcs/carton

**GM062-104B SPES Waterproof Boot cover** SPES 防水靴套

**Specification:**

Color: white, blue, green, etc.  
Material: SPES fabric  
Comfortable elastic hem, helps you keep your cool.  
High water proof efficiency and high repellency  
Weight: 30G/35G/40G/45G.. etc.  
Standard size: 40 x35cm, or other size per customer' demand

GT062-104B



**GM063-100 Spunbond Meltdown Spunbond( SMS ) Surgeon Gown** SMS 手术衣

**Specification:**

Material: Spunbond Meltdown Spunbond( SMS )  
Color: white/blue/yellow/green.  
Weight: 40gsm/ 45gsm/48gsm/50gsm/55gsm etc.  
Size: S/M/L/XL/XXL  
Suited to hospital use & for surgery

**Specification:**

With ties at back of the neck and waist  
Comfortable sleeve tailoring for maximum freedom of movement  
Soft knit cuffs for wrist comfort  
Secure sleeve seams for protection  
Extra protection reinforcement on sleeve & chest and full-body on selected gowns  
Spunbond meltdown spunbond(SMS) material gown, breathable, quiet and supple, high repellency, AAMI standard, sterile with comfortable design for operating theater  
With knitted cuff, Velcro neck back closure, side waist band  
1 hand towel is available  
1 wrap towel is available  
Sterile method: non sterile or EO STERILE,  
Packing: 10 pcs/poly bag or 1 pc/sterile bag



**GM063-101 Reinforced SMS Surgical gown** 加强型 SMS 手术衣

**Specification:**

Material: Spunbond Meltdown Spunbond( SMS )  
Color: white/blue/yellow/green.  
Weight: 40gsm/ 45gsm/48gsm/50gsm/55gsm etc.

Size: S/M/L/XL/XXL

Suited to hospital use & for surgery

**Specification:**

With ties at back of the neck and waist  
Comfortable sleeve tailoring for maximum freedom of movement  
Soft knit cuffs for wrist comfort  
Secure sleeve seams for protection  
Extra protection reinforcement on sleeve & chest and full-body on selected gowns  
Spunbond meltdown spunbond(SMS) material gown, breathable, quiet and supple, high repellency, AAMI standard, sterile with comfortable design for operating theater  
With knitted cuff, Velcro neck back closure, side waist band  
Reinforced arm and additional liner at front  
1 hand towel is available  
1 wrap towel is available  
Sterile method: non sterile or EO STERILE,  
Packing: 10 pcs/poly bag or 1 pc/sterile bag



**GM063-102**  
**SMMS Surgical gown** SMMS 手术衣

Fluid repellent with a double layer of meltblown material to improve performance. Gowns feature an overlap sterile back and adjustable velcro fastening and soft fabric cuffs for extra comfort.

**Specifications:**

- 5-layer anti-static breathable fabric
- Fluid repellent
- Excellent microbial barrier
- Excellent physical barrier to contaminants and fluids
- Stronger fabric with a matte finish, resistance to tears
- Softer fabric with comfort inside and outside
- Generous cut
- Velcro neck binding
- Book fold.

Can be used in hospital, laboratory, food industry, any where need protection.

Prevent and isolate dust, particle, alcohol, blood, bacte

Packed with 2 sterile hand towels

**Specification:**

Material: SMMS

Color: white/blue/yellow/green.

Weight: 40gsm/ 45gsm/48gsm/50gsm/55gsm etc.

Size: S/M/L/XL/XXL

Suited to hospital use & for surgery

Packed with 2 sterile hand towels

Packed with 1 wrap towel

Sterile method : EO STERILE,

Packing: 1 pc/sterile bag



**GM063-110**  
**Disposable Isolation Gown** 一次性PP隔离服

**Specification:**

Material: offer economical, comfortable reliable and protection, Non woven PP Fabric.

Full length for maximum coverage, Latex free , Sewn seams provide superior strength,

and Extra length on waist ties allow gown to be easily secured in front With ties at back of the neck and waist, Design for the best comfort of all medical personnel, cost effective, with great protection

Color: white/blue/yellow/green.

Weight: 18/ 25/30/40GSM or 18/ 25/30/40GSM PP + 15 PE

Size: 115 x137cm or 120x140cm etc.

A. Elastic wrist, B. with knitted cuff

10pcs/bag, 100pcs/ Carton or 1 pc/sterile pack ; 50pcs/carton

Sterile method: EO STERILE

Suited to Lab, industry, hospital and daily home use.

**Form of packing:**

10 pcs/poly bag or 1pc/sterile packing

Carton size: 560x320x350mm



**GM063-103**  
**Spunlace Surgical Gown** 水刺布手术衣

Anti-bacterial, alcohol-repellent, blood-repellent Waterproof dustproof breathable, much soft, comfortable mobility, high repellency, AAMI Standard, Made from soft spunlace

Size: S/M/L/XL/XXL

1 hand towel is available

1 wrap towel is available

Sterile method : non sterile or EO STERILE,

Packing: 10 pcs/poly bag or 1 pc/sterile bag

Packing: 1 pc/sterile bag, 25pcs. carton,

Carton size: 340x280x380mm,



**GM063-120**  
**PP+PE Isolation Gown** 一次性 PP + PE 隔离服

**Specification:**

Material: offer economical, comfortable reliable and protection, Non woven PP+PE Film(100% waterproof efficiency )

Full length for maximum coverage, Latex free ,

Sewn seams provide superior strength,

and Extra length on waist ties allow gown to be easily secured in front With ties at back of the neck and waist

Color: white/blue/yellow/green.

Weight: 18/ 25/30/40GSM or 18/ 25/30/40GSM PP + 15 PE Film

Size: 115 x137cm or 120x140cm etc.

Elastic wrist or with knitted cuff

10pcs/bag, 100pcs/ Carton or 1 pc/sterile pack ; 50pcs/carton

Sterile method: EO STERILE

Suited to Lab, industry, hospital and daily home use.

**Form of packing:**

10 pcs/poly bag or 1pc/sterile packing

Carton size: 560x320x350mm





Heat Sealed Seams

**GM063-130 Impervious Isolation Gown** 一次性超声波隔离服

**Specification:**  
Material: offer economical, comfortable reliable and protection, Non woven PP/ PP+PE  
100% waterproof efficiency and high repellency, AAMI Standard, Heat Sealed Seams by automatic Ultrasound machine, With Ties at back of neck ,  
Extra length on waist ties allow gown to be easily secured in front

Color: white/blue/yellow/green.  
Weight:18/ 25/30/40GSM or 18/ 25/30/40GSM PP + 15 PE  
With Elastic wrist or with knitted cuff  
Size:115 x137cm or 120x140cm etc.  
10pcs/bag, 100pcs/Carton or 1 pc/sterile pack ; 50pcs/carton

**Form of packing:**  
10 pcs/poly bag or 1pc/sterile packing  
Carton size: 560x320x350mm



**GM063-150 CPE gown** CPE 隔离服

This kind of CPE gown Can be used in hospital, laboratory, food industry, any where need protection.


**Descriptions:**  
Can be used in hospital, laboratory, food industry, any where need protection. Prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading. Can be used as a visit coat .Prevent te cross infection in medical treatment at sanitation field. Economical and disposable. Accord with the hygienic standard.

- 1.Material:CPE
  - 2.Characteristics: prevent and isolate dust, particle, alcohol, blood, bacteria and virus invading.
  - 3.Thumb loop wrists prevent sleeves from sliding up when putting gloves on
  - 4.Colour:blue
- Size: 125x200cm, 62g.  
Package: 1 pc/bag, 100pcs/carton  
Carton size: 460x200x250mm



**GM063-141 Disposable Exam Paper Towel** 一次性擦手巾

Made from a high quality soft paper material , towels are unembossed and can be used as disposable premium hand towel for surgery.  
Size: 45cm x45cm, 50cm x60cm or other size per customer's demand



**GM063-142 Disposable Exam Table Paper** 一次性检查纸卷

Our Exam Table Papers provide gentle patient comfort. Durable protection for any size exam table. Choose from sturdy Smooth Paper, or select soft Crepe for a quieter comfort.  
Size: 50cm x 50 mtrs , or Large sizes are perfect for bariatric needs.



**GM167-101 Sterile Delivery Pack** 一次性消毒产包

Non-woven surgical gown-1PC;  
86x80cm Packing Clothes 2ply(Non-woven and thin film)-1PC;  
97x73cm Sheet 3ply(Thin film,non-woven and paper pulp)-1PC;  
97x73cm Triangle Trouser 3ply(Thin film,non-woven and paper pulp)-2PCS;  
66x43cm Non-woven Cloth 2ply(Thin film and non-woven)-2PCS;  
66x43cm Non-woven Cloth 1ply-1PC;Umbilical Cord Clamp-1PC;  
8cm Cotton sticker-2PCS;PE glove-1PAIRS;  
Latex examination glove-1PAIRS;  
8x10cm Gauze swab-5PCS;Face mask-1PC;Cap-1PC



**GM063-301 Lab Coat** 一次性实验服

**Specification:**  
Material: Non Woven PP/PE  
Color: White/Blue/Yellow/Green.  
Weight: 30/35/40/45GSM ,  
Size: M:105x130 cm / L:105x135 cm  
XL:115x140 cm / XXL:115x145 cm  
With Label Collar/ Knitted Collar  
and 4 Press Buttons  
Button Closure or Zipper Closure  
or Velcro Closure  
Size: M/L/XL/XXL  
Suited To Hospital and Daily Home Use. Form of Packing:1 0Pcs/ Inner Poly Bag; 100Boxes/Carton  
Gross.Weight:12.5kgs  
Net.Weight:11.0kgs



**GM063-302 Lab Coat** 一次性实验服

**Specification:**  
Non Woven PP Fabric 40g, 45g  
71x62cm Short Sleeve Top  
4 Button Front Closure  
108x56cm Pant With Elastic Band  
Color: Blue  
Keeping The Environment Free of Infection and Disease Is Always a Battle  
For a Health Professional,  
Suited To Hospital and Daily Home Use. Form of Packing:1Pcs/ Inner Poly Bag ; 50Boxes/Carton  
**Size of Packing:**  
Export Carton: 560x350x500mm  
Gross.Weight: 12.5kgs  
Net.Weight:11.0kgs



**GT063-400**  
**Disposable Scrub Coat**  
一次性洗手衣

**Description:**  
Disposable Scrub Coat are Perfect For any Medical Professional Shopping on a Budget. Available In Several Colors And Made From Quality Non Woven PP Or SMS Material, Disposable Scrub Coat Are Attractive As Well As Affordable. Suited To Hospital And Daily Home Use.

**Specification:**  
Non Woven PP Or SMS Fabric:48g,55g  
V Collar Top  
Elastic Waist Pant  
Color: Green, Blue  
Unisex  
Form Of Packing:1Pcs/ Inner Poly Bag ;  
50Boxes/Carton  
**Size of packing:**  
Export carton:430x340x550mm  
Gross.Weight:8.5kgs  
Net.Weight:7.5kgs



**GT166-100 Body bag** 裹尸袋

10pcs/polybag;  
20bags/ctn  
42x34x40/200pcs



**GT160-100**  
**Non-woven Fabric** 无纺布样

**Specification:**  
PP, SMS, Spunlace  
Kinds of Gram  
Kinds of Color



**GT057 -101**  
**Doctor coat** 医生服

- officer collar
- fly front, full length with 5 button
- 1 left breast pocket and 2 slanted pockets
- center back vent, waist with back belt
- size: 34 " -56 " , could specially make
- 65% polyester , 35% cotton twill or 100% combed cotton fabric
- with white color only



**GT057 -102**  
**Lab coat** 实验服

- peaked lapel collar
- 4 button closure
- 1 left chest pocket with pencil division
- 2 roomy patch pockets
- knee length
- size: 32 " -60 " in regular and long length is available
- 65% polyester , 35% cotton twill or 100% cotton fabric
- color: white, green, blue



**GT063-500 Overall Gown**  
一次性防护服

**Specification:**  
With Head Hood & Shoe Cover and Zipper  
Material: PP/ PP+PE Film / SMS  
Color: White/Blue/Green/Dark  
Weight: 35/45/50/60GSM  
**Size:** M/L / XL/XXL at Clients'option  
S: 164x126 cm / M: 172x130 cm  
L: 178x134 cm / XL: 182x138 cm  
XXL: 184x142 cm

Suited To Lab, Industrial, Agriculture, Hospital And Daily Home Use.  
Unique One Piece Design  
That Cover Whole Body Convenience And Easy To Wear And Off  
Available In Standard or Pe Coated Coverall, Zipper With Storm Flap, Attached Hood, Boots, Elastic Wrists,  
Waist & Ankles Are Available  
Bound Seam Construction  
Packing: 1Pcs/Poly Bag 50Bags/Carton  
Carton size: 480x340x410mm

**GT057 -103**  
**Unisex solid color scrub** 洗手衣

- V-neck pull over scrub top**
- short sleeves
  - 1 left breast pocket and 2 patch pockets and side slits
- Unisex pant**
- elastic waist band
  - slash pockets
  - front fly, zipper closure

**65/35 T/C twill or 100% cotton fabric**  
**size: XS, S, M, L XL, 2XL**  
**color: white, teal, navy, royal blue, khaki, hunter green, purple solid**  
**could be tailored as per any patterns**





**GM057 -105**  
**Nurse skirt 护士服**

- label pocket
- pull over without other closure
- 2 roomy pocket
- Level stitch above the breast
- 100% cotton combed fabric

**GM057 -104**  
**Doctor shirt 医生服**

- Round collar
- Shorter sleeve
- 1 left breast pocket and 2 patch pockets (with 4cm top hem)
- 6 button closure
- 100% cotton or linen poplin
- Color : white, teal
- Size:XS,S,M,L,XL,2XL



**GM057 -106**  
**Nurse skirt 护士服**

- Round collar
- 5 button front closure
- color mends for decorated
- 65/35 T/C multi strip poplin

**GM057 -107**  
**Nurse skirt 护士服**

- Sweetheart collar
- pull over without other closure
- 2 patch pocket
- 100% cotton combed fabric



**GM062-104**  
**High Elastic shoe cover 高弹鞋套**

Made from High Elastic non-woven Fabric.  
Size: H23.1 X W17mm  
Color: White,blue,green,red,brown,yellow Etc.  
Package: 260 pcs/bag, 6 bags/carton  
carton size: 64x47x45 cm



**GM168-103**  
**Disposable string 一次性无纺短裤**

Non Woven PP, 35 g/m2  
one size (fits all).  
Package: 100 pcs/bag, 10 bags/carton  
carton size: 46x30x29 cm



**GM168-104**  
**Disposable pant 一次性无纺短裤**

Non Woven PP, 35 g/m2  
one size (fits all).  
Package: 100 pcs/bag, 10 bags/carton  
carton size: 46x30x33 cm




**GM079-100 O. R. Towel 手术巾**

Material:100% Cotton,  
Pre-washed or Non-washed.  
Standard Weight: 200Gram/m<sup>2</sup> or 250Gram/ M<sup>2</sup>  
or Other Gram Per Customer' Demand.  
Color: Blue  
Stardard Opened Size: 43x66cm(±2cm) or Other  
Size Per Customer' Demand  
Packing: Sterile Package or In Bulk  
10Pcs / Poly Bag, 20Bag / Ctn,  
Dim: 46 X 44 X 33cm, 15.5/15Kgs



**GM078-100 洗手刷**  
**Surgical Hand Brush**  
with / without Iodine



**GM064-100 Drape With Hole 消毒手术洞巾**

Material: PU or Lo Drape  
Pu Drape Is a Transparent, Adhesive Polyurethane Film,  
Which Can Adhere Throughout Long Surgical Operations To The Surrounding Skin And  
Most Importantly To The Wound Edge, Thereby Significantly Lessening The Chance Of  
Compromising The Sterile Environment.  
Iodrape Is An Enhanced Surgical Drape Product That Is Impregnated With Iodophor P.V.P.  
Iodrape Helps To Maintain A Sterile Barrier, Lowers The Risk of Microbial Transmission,  
and Significantly Lowers The Colony Forming Units (CFU'S) of The Bacteria Staphylococci  
During Use In Operations. The Need For The Enhanced Product Is Especially Useful For The  
Following Operations: Orthopedic Procedures, Cardiac Valve Replacements, Venous  
Retention Suture <needle>, Peritoneal <encheiresis> Therapy, And By-pass Operations  
of The Cerebrospinal Fluid.

**Features:**

- Breathable
- Transparent and Thin
- Conformable and Elastic
- Easily Applied To The Skin
- Sterilization: EO
- Size: 50cmx50cm, 70cmx90cm,70cmx110cm,  
150cmx200cm or Other Size Per Customer' demand




**GM168-101 Underpad 一次性床垫**

Material: Non-Woven fabrics + FULLF +SAP + PE Film  
Description: Color: blue  
60X60 cm 10pcs/bag 10bags/ctn 57X22X31 cm  
60X90 cm 10pcs/bag 10bags/ctn 42X28X42 cm



**GM064-101 Non-Woven fabrics Bedsheet 一次性无纺床单**

**Description:**  
Soft touch and harmless to your skin  
100% biodegradable and environmentally friendly  
SPEC: light blue, 60x180cm,  
30gsm other size and specification are also acceptable  
PACKING: 25pcs/polybag; 4bags/ctn  
CTN SIZE : 32x25x24 cm /100pcs



**GM064-102 Non-Woven fabrics Bedcover 无纺床罩**

**Description:**  
disposable bed cover made of pp nonwoven fabric,  
helps you to be clean & hygiene on your travel and regular purpose.  
SPEC.: 30g,210x90x20cm,  
other size and specification are also acceptable  
PACKING: 1pc/roll;10rolls/bag, 2bag/ctn  
CTN SIZE: 31x14x26cm /20pcs



**GM168-102 Adult Diaper 成人纸尿裤**

**Specification:**

65x80cm, with leg cuffs, waistline:71-121cm  
65x80cm, without leg cuffs, waistline:71-111cm  
80x96cm, with leg cuffs, waistline:96-147cm  
80x96cm, without leg cuffs, waistline:96-137cm

**Material:**

1. Inner layer: Polyolefin Non-Woven
2. Absorption :: Cotton Pulp, Tissue, SAP
3. Stopper :: Adhesive tape, Frontal tape
4. Back sheet: Water-proof PE Back Layer
5. Other : Leak guard, Fabric Glue and so on

**Features:**

- 1) Wetness indicator
  - 2) Super absorbent core, dry and comfortable
  - 3) Refastening Frontal tape is convenient to put on and removed
  - 4) PE film makes the surface dry and comfort
  - 5) Super soft non-woven surface
  - 6) Leak Guard prevent leakage efficiently
- Package: 10pcs/bag, 80pc/carton  
Carton size: 680x255x480mm



**GM064-103 Nonwoven Pillow cover 无纺枕套**

Disposable pillow covers are used widely in hotels, hospitals, home, and airline and in all travels


**Specification:**

80gsm,50x60cm,with 5cm gusset,  
other size and specification are also acceptable  
PACKING: 100pcs/polybag; 1bag/ctn  
CTN SIZE: 54x32x38 /100pcs



DISPOSABLE  
SURGERY  
SET





**GMCP-001 Universal Set 通用包**

- 1 Table Cover 140 x 190cm
- 1 Mayo Stand Cover 80 x 145cm
- 1 Bottom Drape 160 x 180cm
- 2 Side Drape 75 x 98cm
- 1 Top Drape 160 x 240cm
- 4 Hand Towels 30x 40cm
- 2 Op Tape 10 x 50cm



**GMCP-005 膝关节手术包 Stand Orthopedic Drape Packs for Extremity Plastic Surgery**

- ORTHOPEDIC SURGICAL DRAPE PACK
- 1 Table cover 150 x 300cm
  - 1 Mayo stand cover 63 x 152cm
  - 1 Imervious split drape 150 x 190cm
  - 1 Lower extremity drape 260 x 380cm
  - 1 Untility drape 50 x 75cm
  - 2 Medium drape 150 x 200cm



**GMCP-008 Caesarian Delivery Pack 产包**


- 1 Back Table Cover, Reinforced 140 x 190cm
- 2 Absorbent Towel 40 x 65cm
- 1 Mayo Stand Cover, Reinforced 75 x 145cm
- 4 Utility Drapes with adhesive tape 75 x 50cm
- 1 Utility Bag
- 1 Whole body drape T shaped 255 x 300 x 190cm fenestrated with control fabric reinforcement 30 x 30cm with tube holder on both holders with fluid collection pouch with suction por L/S
- 2 umbilical cord clamp
- 1 needle counter
- 2 adhesive tape (steri-tape) 10 x 50cm
- 1 Baby blanket 80 x 75cm
- 1 placeta Basin
- 1 bulb syring (ear ulcar)

**GMCP-009 Dental Pack 牙科包**

- 1 support for tubing with adhesive
- 1 fenestrated drape 75 x 75cm
- 1 surgical drape 100cm x 100cm
- 1 surgical drape 150 x 150cm
- 2 tubing protection 120 x 80cm
- 2 isolation gowns
- 2 scialytic lamp handle cover
- 2 surgical face mask with ties
- 1 nurse cap with elastic
- 2 surgeon cap with ties
- 1 nail brush
- 2 hand towel
- 4 polyethylene glove
- 1 saliva ejector

**GMCP-1523 Laparotomy Drape Pack 腹腔镜检查手术包**

- 1. Laparoscopy Drape 250 x 300cm Main drape size: 150 x 250cm, B190 x 203cm x 1pc with adhesive reinforcement size: 73 x 73cm
- 2. Large Table Cover 150 x 200cm x 1pc
- 3. Small Table Cover 80 x 150cm x 1pc
- 4. Op-tape 10 x 30cm x 1pc
- 5. Surgical Gowns reinforced x 2 pcs
- 6. Hand Towels 40 x 40cm x 3 pcs
- 7. Universal Bag x 1pc
- 8. Sterile Marker x 1pc








**GM170-100 Knee Arthroscopy Set 膝关节镜检查手术台组件**


- 1 pc bilplex drape mm. 2000x2000;
- 1 pc mayo cover mm. 800x1400;
- 2 pcs towel mm. 400x400;
- 2 pcs Adhesive strip mm. 100x500;
- 1 pc orthopaedic stockinette mm. 220x750;
- 1 pc knee arthroscopy drape (Double elastic aperture, Fluid collection pouch with suction port);
- 1 pc Triplex Drape mm. 1500x2000 with adhesive edge.

**GM170-200 Ophthalmology drape 眼科手术用铺巾**

Base cloth of spun laced nonwoven faric is adopt nonwoven faric cloth which contains less hemp(low dust).








**GM170-700 Cardiovascular Set**

心血管手术台组件

1pc bplex drape mm.2000x2400;  
1pc mayo cover mm.800x1400;  
4pcs towel mm.400x400;  
4pcs Adhesive strip mm.100x500;  
2pcs bplex drape mm.750x1000 with adhesive edge;  
1pc cardiovascular drape (full chest and femoral incise area. Adjustable tube holder. Instrument pouches.)

**GM170-800 Cardiac surgical drape** 心脏外科手术用辅巾

Water-repellent spun laced nonwoven fabric which is the same as base cloth shall be used as large shielding material and installed after intergrating with the base cloth. In addition, surrounding parts of survey section shall be reinforced with water-absorbable/water proof material.



**SMS DRAPES** SMS 手术辅单

<b>GM1801</b> small drape 40cmx60cm	<b>GM1802</b> half drape 98cmx152cm	<b>GM1803</b> medium drape 112cmx190cm	<b>GM1804</b> large drape 150cmx190cm
<b>GM1805</b> 60x60cm 5x9cm egg shaped fenestration	<b>GM1806</b> 80x80cm 6x17cm fenestration	<b>GM1807</b> 94x94cm	<b>GM1808</b> 150x80cm 5x15cm Fenestration
<b>GM1809</b> 152x110cm 10x15cm Fenestration	<b>GM1810R</b> 120x150cm 30x35cm reinforced 10x15cm Fenestration	<b>GM1811</b> 124x180cm 10x15cm Fenestration	<b>GM1812</b> 137x183cm 7x16cm Fenestration
<b>GM1813</b> 170x150cm 5x15cm Fenestration	<b>GM1814</b> 300x100cm 20x15cm Fenestration		

**OPHTHALMIC DRAPES** 眼科手术单

<b>GM0301A</b> 40x60cm fluid pouch(+support strip) 4.4x7.2cm incise Fenestration	<b>GM0301B</b> 40x60cm fabric+PE W/ Fenestration	<b>GM0302</b> 101x127cm fluid pouch(+support) 9x11cm incise
<b>GM0303A</b> 140x140cm fluid pouch(+support strip) 9x11cm incise	<b>GM0303B</b> As per PM-0304 +4.4x6.5cm oval incise	<b>GM0304</b> 158x168cm dual fluid pouches 9x12cm incise
<b>GM0305</b> 200x140cm fluid pouch(+support strip) 9x11cm incise	<b>GM0306</b> 180x250cm dual fluid pouches(+3support strips) 9x12cm incise with 4.4x6.5 oval Fenestration	

**LITHOTOMY DRAPES** 结石手术单

	<p><b>GM1201</b> Tdrape 150x190cm 60x26cm reinforcement 8x12cm Fenestration</p>		<p><b>GM1202</b> T drape large 150x220cm 60x26cm reinforced 8x12cm Fenestration</p>		<p><b>GM1203</b> Tdrape with mesh 150x220cm 60x26cm reinforcement 8x12cm Fenestration</p>
--	---	--	---	--	---

**EXTREMITY DRAPES** 四肢手术单

	<p><b>GM1401A</b> Unilateral drape 210 x 320cm 80 x 115cm reinforcement 6 1/2 cm expandable Fenestrations cord holders (as shown)</p>		<p><b>GM1401B</b> Bi-lateral drape 210 x 300cm 80 x 110cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown)</p>
	<p><b>GM1402</b> Unilateral drape large 260 x 300cm 80 x 110cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown) Side view panels</p>		<p><b>GM1405</b> Hand drape 300 x 295cm 40cm x 86cm reinforcement 6 1/2 cm expandable Fenestrations</p>
	<p><b>GM1403A</b> Unilateral drape 224 x 320cm 80 x 110cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown) Side view panels</p>		<p><b>GM1403B</b> Bi-lateral drape 240 x 300cm 80 x 100cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown) Side view panels</p>
	<p><b>GM1404A</b> Unilateral drape+armboards 320 x 290cm 80 x 110cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown)</p>		<p><b>GM1404B</b> Bi-lateral drape+armboards 320 x 290cm 110 x 120cm reinforcement 6 1/2 cm expandable Fenestrations Cord holders(as shown)</p>

**DENTAL DRAPES** 牙科手术单

	<p><b>GM0501</b> 150 x 195cm 96 x 74cm reinforcement 10 x 10cm fenestration</p>		<p><b>GM0502</b> 180 x 250cm 70 x 70cm reinforcement 10 x 10cm fen fenestration</p>		<p><b>GM0503</b> 180 x 270cm 4 x 13cm fenestration</p>
--	---	--	---	--	--

**HICKMAN DRAPES** 颈部手术单

	<p><b>GM0601A</b> Hickman with oval 170 x 270cm 80cm x90cm reinforcement 15x 20cm oval fenestration with d/s tape</p>		<p><b>GM0601B</b> Hickman with arch 170 x 270cm 80cm x 90cm reinforcement 17x 20cm arch fenestration with d/s tape</p>
--	---	--	--

**PROCEDURAL DRAPES** 手术单

	<p><b>GM2101A</b> 60 x 60cm 15cm x 15cm d/s tape over 6cm fenestration</p>		<p><b>GM2101B</b> 60 x 60cm 15cm x 15cm d/s tape over 6cm fenestration Perforation edge to centre</p>
--	--	--	---

**BRACHIAL DRAPES** 手术单

	<p><b>GM1501TR</b> BRACHIAL DRAPE 72 x 100cm 7.5 x 91/2cm fenestration over d/s tape.</p>		<p><b>GM1502</b> 90 x 120cm 30 x 40cm reinforcement 6 1/2cm incise over 12 1/2cm fenestration</p>
--	---	--	---

**HEAD DRAPES** 头部手术单

	<p><b>GM0101</b> Head drape 140 x 86cm L/blue 114 x 64cm D/blue</p>
--	---

**MISCELLANEOUS DRAPES** 多用手术单

	<p><b>GM2201</b> 90 x 100cm 30 x 40cm Deep arch</p>		<p><b>GM2202</b> 120 x 150cm 10 x 30cm fenestration</p>
--	---	--	---