



# GlenMedicon LLC Medical Instruments Co.,Ltd.

*Best Greetings to All Patients in the World!*

▶ Medical PVC Tube Version 2019

[www.glenmedicon.com](http://www.glenmedicon.com)



**GLENMED**  
*- we help cure -*

## GlenMedicon LLC

Add: 7122 S Sheridan Rd Ste 2 - 1012Tulsa, OK 74133-2748

Tel: 1-918-8091339

Fax: 1-918-523-7121

E-mail: [info@glenmedicon.com](mailto:info@glenmedicon.com)

Http://[www.glenmedicon.com](http://www.glenmedicon.com)



[www.china-greetmed.com](http://www.china-greetmed.com)



**GlenMedicon** Medical Equipment & Supplies is a minority and women owned business enterprise. We provide practitioners, patients and the general public with the highest quality products at the lowest available prices. We envision a world in which medical products can be accessible to all.

Regardless of the market or product category, our diversified relationships with all CE approved and FDA approved product manufacturing partners around the globe, can help your company drive revenue and overall valuation. Our strong track record of proven performance is the direct result of our decision to build the company's future on a specialized approach to marketing - one that promotes a unique understanding of the needs of the medical community. Our global accessibility gives us resources to effectively fulfill our commitment to excellence, now and in the years to come.

**GlenMedicon** successfully represents some of the most well respected Speciality Medical Device manufactures in the world today, allowing us to expand our product offering to multiple departments in the Hospitals, Diagnostic Laboratories , Rehabilitation Centers, Home health & Assisted living Centers as well as to Surgery Centers and a variety of Specialist Clinics. We specialize in respiratory, urology, wound care , catheter , syringes and needles.

Our primary concern is your satisfaction.

Best Greetings to all the patients in the world.

We pray no disaster ,no flu, We hope everyone is healthy and fortunate.



## **GM017-100A**

### **Suction Catheter (Plain Type)**

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes
- (4) Connector(color coded)-Size(Fr/Ch):  
F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



### **GM017-100B**

#### **Suction Catheter - Thumb Control ( Cap-cone Type )**

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes
- (4) Connector(color coded)-Size(Fr/Ch):  
F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



### **GM017-100C**

#### **Suction Catheter (T Type)**

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch):  
F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.





**GM017-100D**  
**Suction Catheter (Y Type)**

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Color: Transparent or translucent(surface fog form).
- (3) Distal end is opened with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch):  
F5, F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.

**Options Available:**

**More connector:** Transparent connector

**More tip:** declined tip, big tip and so on.

With/without X-ray detectable thread.

**Length:** 45cm-55cms or any length meetcustomers' demands.

Or DEHP free tube.



**GM017-200**  
**Extractor Mucus**

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free available.
- (2) Suction tube length: 35-40cms.
- (3) Capacity: 25ml.
- (4) Size: Fr8, Fr10, Fr12, Fr14, also available as per customer request
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



### **GM023-200**

#### **Closed Suction Catheter**

- (1) Made of non-toxic, non-irritant PVC(Medical grade).
- (2) Size(Adult and pediatric): 6Fr, 8Fr, 10Fr, 12Fr, 14Fr, 16Fr.
- (3) Package: Single blister package.
- (4) Sterilized by EO GAS.



### **GM018-100**

#### **Feeding Tube**


- (1) Made of non-toxic, non-irritant PVC(Medical grade) , or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent(surface fog form)
- (3) Closed tip with two lateral eyes
- (4) With/without X-ray detectable thread
- (5) Connector(color coded)-Size(Fr/Ch): F4,F5, F6, F8, F10, F12, F14, F16, F18
- (6) Two types connectors available.
- (7) Package:Single in PE or blister package
- (8) Sterilized by EO GAS



### **GM019-100A**

#### **PVC Stomach Tube**

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free.
- (2) Color: Transparent or translucent (surface fog form)
- (3) With/without X-ray detectable thread.
- (4) Connector(color coded)-Size(Fr/Ch):  
F6, F8, F10, F12, F14, F16, F18, F20,F22.
- (5) Package:Single in PE or blister package.
- (6) Sterilized by EO GAS.



### **GM019-100B**

#### **100% Silicone Stomach Tube**

- (1) Made of non-toxic, non-irritant 100% silicone.
- (2) Color: Transparent or translucent (surface fog form)
- (3) With/without X-ray detectable thread.
- (4) Connector(color coded)-Size(Fr/Ch):  
F6, F8, F10, F12, F14, F16, F18, F20,F22.
- (5) Package:Single in PE or blister package.
- (6) Sterilized by EO GAS.



**GM019-100C**  
**Duodenal Tube**

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent(surface fog form).
- (3) With X-ray detectable thread.
- (4) Closed distal end with four lateral eyes.
- (5) Connector(color coded)-Size(Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (6) Package: Single in PE or blister package.
- (7) Sterilized by EO GAS.



**GM019-100D**  
**Ryle's Tube**

- (1) Made of non-toxic, non-irritant PVC(Medical grade), or DEHP free or 100% silicone.
- (2) Color: Transparent or translucent(surface fog form).
- (3) With X-ray detectable thread.
- (4) With metal ball.
- (5) Connector(color coded)-Size(Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22.
- (6) Package: Single in PE or blister package.
- (7) Sterilized by EO GAS.



### **GM020-100**

#### **Rectal Tube**

- (1) Made of non-toxic, non-irritant PVC (Medical grade) or DEHP free.
- (2) Color: Transparent or translucent. (Surface fog form).
- (3) Closed tip with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch): F20, F22, F24, F26, F28, F30, F32, F34, F36.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



### **GM021-100**

#### **Nelaton Catheter**

- (1) Made of non-toxic, non-irritant PVC (Medical grade) or DEHP free.
- (2) Color: Transparent or translucent. (Surface fog form).
- (3) Closed tip with two lateral eyes.
- (4) Connector(color coded)-Size(Fr/Ch): F6, F8, F10, F12, F14, F16, F18, F20, F22, F24.
- (6) Nelaton or tiemann, male or female available.
- (7) Package: Single in PE or blister package.
- (8) Sterilized by EO GAS

